

Flavour & Fragrance Natural & Synthetic Compounds GC/MS Library

# FFNSC3 Library

GC/MS is commonly used for flavor and fragrance analyses, and searching a library of GC/MS spectra is convenient for identifying unknown components. However, similarities in the mass spectra of some fragrant compounds makes identification difficult. The FFNSC3 library contains retention indexes using three types of columns as well as mass spectra. Shimadzu's GCMSsolution library search, utilizing the retention index, provides highly accurate identification results and makes it ideally suited for flavor and fragrance analysis.



GCMS-QP2020

## 3462 Mass Spectra

The FFNSC 3 library is registered with 3462 natural and synthetic chemical compounds relating to flavors, providing even greater support for analysis.

## Retention Index Compatible with Three Types of Columns

The FFNSC 3 library contains retention indexes compatible with three types of columns including wax columns. Shimadzu's GCMSsolution library search, utilizing the retention index, narrows down the search results and makes it easy to identify various fragrant compounds.

## FFNSC 3 Library

### Specifications

Registered Compounds:

3462 flavor and fragrance compounds

Registered Retention Indexes:

Micropolar column (SLB™-5ms) 3462[n-alkane], 2516[FAMEs]

Highly-polar column (SUPELCO WAX™ 10) 1466[FAMEs], 1466[FAEEs]

Non-polar column (Equity™-1) 646[n-alkane]

Note 1: Items in parentheses ( ) are Recommended Columns, items in brackets [ ] are reference compounds.

Note 2: SLB, SUPELCO WAX and Equity are trademarks of Sigma-Aldrich Co. LLC

Registered Information:

Mass spectrum, Retention index for each column, CAS number, compound name, molecular weight, compositional formula

Applicable models: GCMS-QP Series + GCMSsolution Ver.2.6 or later

GCMS-TQ Series + GCMSsolution Ver.4.0 or later

(single-quad mode only; does not include MS/MS spectra)

## High Hit-ratio Mass Spectra

The FFNSC 3 library was developed using the Shimadzu GCMS-QP series by the skills-based group led by Prof. Mondello, (Messina Univ. in Italy). Performing analysis using Shimadzu GCMS enables obtaining mass spectra with similarities close to those registered in the library.



The mass spectra in the FFNSC library were obtained by the group led by Prof. Mondello, University of Messina, Italy. The FFNSC library copyright is owned by Chromaleont S.r.l., a spin-off company of the University of Messina.

# FFNSC3 Library

## Filtering with retention index

Since flavor and fragrance components sometimes contain compounds having similar chemical structures, multiple compounds with similar mass spectra are listed as candidates when performing a library search using only the mass spectrum. Filtering with the retention index sorts the candidates by retention index, thereby providing highly accurate identification results.

Results of similarity search in mass spectrum

Hit	Similar	Regi	Compound Name	Mol Wt	Formula	Library No
1	96	<input checked="" type="checkbox"/>	Terpinolene $\text{C}_{10}\text{H}_{16}$ Cyclohexene, 1-methyl-4-(	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
2	95	<input type="checkbox"/>	Terpinene <alpha> $\text{C}_{10}\text{H}_{16}$ 1,3-Cyclohexadiene	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
3	90	<input type="checkbox"/>	Carene <delta 2 > $\text{C}_{10}\text{H}_{16}$ Dicyclo[4.1.0]hept 2	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
4	91	<input type="checkbox"/>	Carene <delta-3 > $\text{C}_{10}\text{H}_{16}$ Bicyclo[4.1.0]hept-3	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
5	90	<input type="checkbox"/>	Terpinene <gamma> $\text{C}_{10}\text{H}_{16}$ 1,4-Cyclohexadie	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
6	89	<input type="checkbox"/>	2,2-Dimethyl-5-methylene norbornane $\text{C}_{10}\text{H}_{16}$	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
7	88	<input type="checkbox"/>	Sylvestrene <iso-> $\text{C}_{10}\text{H}_{16}$ Cyclohexene, 1-met	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib
8	88	<input type="checkbox"/>	Ocimene <(E)-, beta-> $\text{C}_{10}\text{H}_{16}$ 1,3,6-Octatriene,	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib




Results sorted using retention index filtering

Hit	Similar	Regi	Ret. Ind	Compound Name	Mol Wt	Formula	Library No
1	96	<input checked="" type="checkbox"/>	1006	Terpinolene $\text{C}_{10}\text{H}_{16}$ Cyclohexene, 1-methyl-4-(	136	C <sub>10</sub> H <sub>16</sub>	FFNSC 3Lib

## Multi-column compatible

All of the column retention indexes can be confirmed at the same time, which assists in close investigation of identification results.



## Provides greater support for flavor and fragrance analysis

Gas Chromatograph Mass Spectrometer

### GCMS-QP2020

Enables high-sensitivity and high-speed analysis

Single Quadrupole Mass Spectrometer



## Various System Configurations\*

### Headspace Analysis System



### Thermal Desorption Analysis System\*\*



### GC/MS Off-Flavor Analysis System



\* The accuracy of the retention indexes differs depending on the sample introduction method.

\*\* Not available in the U.S.



Shimadzu Corporation

www.shimadzu.com/an/

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.